

Beds and Cabinets

hospital • medical • ward • theatre • specialist • anaesthetic • emergency crash • ICU • medicine • records • utility • dressing • instrument • end of bed single & multi-drawer • linen • bedside cabinets • over-bed tables • wardrobes, below is a small selection of the range

Emergency Trolleys

Overall Dimensions:

Height 900mm
Width & depth 530mm
Drawer sizes 400 W x 445mm D
Pullout work top 450 x 450mm

IV Pole. Adjustment 1750 to 2460mm

5353E.

Top drawer 75mm high
Lower drawer 125mm high

5353E5.

Top two drawers 75mm high
Lower three 125mm high



The **WDROBE** is 430 W x 540 x 1575 mm high

Features:

- One shelf
- Clothes hanger rail with two height positions



Anaesthetic Trolleys

Overall Dimensions:

Height 900mm
Width & depth 530mm
Drawers sizes 400 W X 445mm D
Pullout work top 450 x 450mm
Drawers are fitted with fixed divisions

5353 AT

Top drawer 75mm high
Lower three 125mm high

5353AT5

Top two drawers 75mm high
Lower three 125mm high



Bedside Cabinets

All **DRCAB** models are 440 W x 400 D x 840 mm high
Drawers are 360 W x 340 D fitted with roller slides
Fitted with four castors

Features:

- 1DRCAB 125 mm drawer a cupboard with one shelf
- 2DRCAB two 125 mm drawers and cupboard
- 4DRCAB two 125 mm and two 175 mm drawers



Medicine trolley

Has an aluminium roller door with a automatic snap lock that can also lock the drawer and the pullout work top. The key is required to unlock only. Push handles on both ends and a drop in lower shelf with a 30mm lip and guardrails. The inserts are removeable for future changes or other dispensing system. The drawer and pullout work top are on double-acting slides.

Overall Dimensions:

Overall height 1370mm
Floor to work top 870mm
Width 710mm
Depth 555mm

6348DT

Has three replaceable pantry type shelves for bottles or packets.

6348DTA

Has three replaceable sliding trays that hold 10, 160 L x 90 W x 65mm D cardboard boxes



Stainless Steel Products—Cambridge Sheetmetals



Dressing Trolley



Flat Deck Trolley



Instrument Trolley



Multiple Drawer Trolley



Plaster Trolley



Cylinder Trolley



Cylinder Cart



Single Shelf Trolley



Resus/Emergency Trolley



Utensil Trolley



Ward Trolley



Bowl Trolley

New Zealand Manufactured—Superb Lasting Quality

Trolleys available in different sizes and configurations

Stainless Steel Products



Apron Rack



Autoclave Tray



Bedpan/Urinal Rack



Cooling Rack



Cylinder Hanger



IV Stand



Footstep Nest



Drug Safe



Supplies Trolley



Soiled Linen Trolley



Meal Trolley Trolley A/B

