

## The Dr Lee range of Electrocardiographs



### Dr Lee ECG-120B

We should like to introduce the Dr Lee range of ECG machines. The ECG-120B is a 12 lead, single channel automatic unit which is very competitively priced and comes fully loaded with features including: 2 line LCD Screen, digital filter to eliminate artefacts, simple one touch operation, an easy paper loading mechanism and built in battery recharger for the optional rechargeable NiCad cell. The printer mechanism is a thermal array type and takes standard ECG paper.


**Optional Accessories:**

- Rechargeable NiCad battery
- Mobile cart
- Dust cover
- Carrying bag
- Patient cable hanger

**Included accessories:**

- Patient cable
- Limb electrodes
- Chest electrodes
- ECG Cream
- Ground lead

**OVER 80 SOLD  
IN THE NZ/PACIFIC REGION!!**



### Dr Lee ECG-310A

The Dr Lee ECG-310A is a top-of-the-line Interpretive ECG machine with all the features and accessories you would expect in a machine over double the price. It has 6/3 channel 12 lead operation and is stacked with powerful software based on the micro-computer analysis including paediatric and post-stress ECG against resting ECG. A concise report of detailed measurement and interpretation acts as quick reference for the physicians in-depth analysis. To assure the accuracy of interpretation and clear ECG printing, digital artifact filters eliminate AC power interference, muscle tremor and baseline wandering. Thermal array printer technology provides a clear and sharp ECG traces and alphanumeric patient data to be read easily.

A comprehensive ECG complex is printed with automatic optimum gain and baseline adjustment at the touch of a button. The 310A check and indicates the paper supply, printer open, improper electrode contact and low battery level. A numeric keypad and on-screen instructions make it easy to enter patient I.D. such as, age, sex, height, weight and blood pressure.

**Included accessories:**

- Patient cable
- Limb electrodes
- Chest electrodes
- ECG Cream

**Optional accessories:**

- Rechargeable NiCad battery
- Mobile cart
- Dust cover
- Carrying bag