

Hadeco vascular Doppler's



Echo Sounder ES-101EX (8MHz)

- Cell-phone size Doppler with built-in speaker & transducer.
- Large LCD for numerical data, PEAK & MEAN velocities.
- 2 AA alkaline batteries.



Minidop ES-100VX

- Pocket Doppler with high quality sound
- Powered by 1 X 9V alkaline battery



Bi-directional Doppler ES-100V3

- Foetal Doppler with large LCD screen for heart rate and monitoring waveform
- Upto 33 minutes of waveform can be stored for each monitoring.
- 30 waveform memory and USB interface
- Powered by 1 X 9V alkaline battery
- Vascular Link Software Smart-V-Link for Windows (Option).
- PPG probe also options



Smartdop® 45

- Vascular Doppler with LCD screen and fast printer
- Upto 33 minutes of waveform can be stored for each monitoring.
- 30 waveform memory and USB interface
- Ni-MH rechargeable battery or AC adapter
- PPG & Pnemo probes (Option).



Echo Sounder ES-102EX

- Desktop Doppler with large speaker
- Displays velocity, heart rate and motion
- Ni-Cd battery pack or AC adapter
- USB optional computer interface.
- Vascular Link Software Smart-V-Link for Windows (Option).
- PPG & Pnemo probes (Option).