

Infusion pumps from GRASEBY & BAXTER

GRASEBY 3300

The GRASEBY 3300 is fully configurable for the most demanding clinical needs and flexible enough for the most demanding patient needs.

Safe The locking syringe cover provides additional security to the standard syringe tampering alarms.

Effective With features like *clinician over-ride*, *loading dose* and a *comprehensive alarm system*, the GRASEBY 3300 is suitable for all types of patient in any clinical situation.



Flexible The 3300 PCA Pump is suitable for use with strong and weak analgesics for both adults and children. A range of syringe brands may be used and a variety of external devices can be connected via the RS232 port.

Easy-to-use The 4-step programme sequence allows for quick and simple programming. The clear LCD display gives the user comprehensive messages and prompts.

Unique The unique Patient Activated Timed Infusion Mode makes a background infusion of morphine safe. This helps to control a patients pain at night to enable better rest recovery.

Accumate 2300

Flow range, ml/hr: 0.1-999

Accuracy, % 5%

Drug/dose calc: yes

Vtbi selector, ml

Total volume infused

Display: ml

Iv set:

Compatibility: open system

Free-flow protection

Air-trapping capability

Needleless iv connection

Front panel lockout

Audible alarm

Volume control

Momentary silence

Nurse call jack

Fluid resistant

Power source, VAC 220volt – 60 hz.

Battery: rechargeable built in

Alarms & indicators

Occlusion upstream

Occlusion downstream

Detection mechanism - selectable

Pressure, psi specify

Real-time display

Infusion complete

Air in line

Empty fluid counter

Flow error

Drop sensor

Set disengaged

Door open

User prompts

Circuit malfunction

Low battery

