



## Water Purification

Crystal Clear Water Filter



Model 724 Water Distiller



The FANEM 724 Water Distiller Size 1 produces 2 litres per hour and the Size 2 produces 5 litres per hour, ideally used in clinics and chemical laboratories.

The REXALL Crystal Clear Water Filter is the ideal companion for your Laboratory, Kitchen, Hospital or Ice Plant. The filter gives the added protection from such water born bacterial diseases including cholera, dysentery, gastritis and typhoid. This economical filter also removes trace elements making it ideal for your film processor, autoclave or steriliser.

Re-generation is simple, just add salt water (1kg salt to 8 litres of water) by removing the top cap, pouring in mixture and then drain.