

Monitoring



The Corometrics 250 Series combines outstanding clinical capabilities with the reliability you demand. These advanced maternal/fetal monitors comprise a high-tech system that allows you to deliver a fulfilling experience for every patient – including the opportunity to celebrate the birth with a choice of 3 songs. Only the best for baby. For every patient, nothing short of excellence will do.

- *The Corometrics 250 Series fetal monitoring capabilities include Fetal ECG for continuous, direct fetal heart rate*
- *(FHR) data with waveform display and printing.*
- *Fetal heart rate high/low alarms provide both audible and visual alerts if the FHR is outside of the user-defined high/low range.*
- *Optional and unique Spectra Alerts™ simultaneously analyze FHR and uterine activity (UA) information, and alert clinicians of deviations from the norm.*
- *Our optional and exclusive keyboard also facilitates basic documentation on the paper strip chart for monitors not connected to a Clinical Information System.*
- *Keeping tabs on mom Maternal well-being is frequently assessed during the physically demanding process of labor, delivery and recovery.*
- *The Corometrics 250 Series provides a choice of integrated GE Ohmeda™, Masimo SET® or Nellcor® maternal pulse oximetry (M_{SpO}2), making it easy for the clinician to monitor, document, and even view the waveform.*
- *The selectable 3-lead maternal ECG with pacemaker detection and rejection generates maternal heart rate, QRS waveform display, and a six-second snapshot printout when requested.*