

## The SOEHNLE range of scales

### SOEHNLE 8310.00 Baby Scale



- Simple one-key operation and tare
- Blankets can be tared
- Removable baby weighing tray: can be used as personal scale up to 20Kg
- Quattro-Tronic: fully scale with 4 load cells for repaid precise weight measurement
- No wear parts and maintenance free
- Easy to read LCD with 20mm high digits
- 9V battery operation with auto-power off power saver.
- Optional power supply
- Kg or lb operation

**Weight Range: 20Kg**  
**Divisions: 10g**

### SOEHNLE 7703 Bariatric Scale



- Approved to class III, MPG Class 1
- BMI, hold size: 600X600mm
- Footplate: 55mm high
- 2 integrated rollers make the scale mobile and easy to move
- Display case made of aluminium extrusion, 180° rotation
- 7 digit LCD backlit 20mm numerals

#### OPTIONS

- Charge Adapter with rechargeable batteries
- Height rod
- Electronic body height meter

**Weight Range: 300Kg**  
**Divisions: 100g**

### SOEHNLE 7800 Adult Digital



- Simple 1-key operation
- Quattro-Tronic with 4 load cells.
- Weighing platform: 310 X 275mm or 300 x 300
- 7 digit LCD backlit 20mm numerals
- Kg or lb

#### OPTIONS

- Charge Adapter with rechargeable batteries
- Height rod
- Electronic body height meter

**Weight Range: 200Kg**  
**Divisions: 100g**

### SOEHNLE 6163 Mechanical scale



- Large, easy to read dial
- Sturdy steel housing
- Classic design
- Weighing platform with easy to clean plastic cover

**Weight Range: 150Kg**  
**Divisions: 1Kg**